



Human B7-H5 / Gi24 / VISTA

Product Name	Human B7-H5 / Gi24 / VISTA
Catalog No	IT-300-043p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged Human B7-H5 / Gi24 / VISTA (a.a.33-194) (GenBank#AAH20568).
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Human VISTA

Note: