

Anti- H5 (PAN)



Product Name	Anti- H5
Catalog No	IT-003-005M28
Source	Mouse monoclonal antibody, clone 2 E10
Isotype	IgG1
Formulation	Each vial contains 100µg purified IgG in 100µl PBS with 40% glycerol
Immunogen	DNA vaccine expressing HA1(A/chicken/Ghana/20/2015(H5N1)) (IT-003-00540p).
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with some H5 as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

