


# Human Poliovirus Receptor-Related Protein 2/PVRL2/CD112



Product Name	Human Poliovirus Receptor-Related Protein 2/PVRL2/CD112
Catalog No	IT-300-078p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged Human Poliovirus Receptor-Related Protein 2/PVRL2/CD112 (aa 32–360; GenBank# Q92692).
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Human Poliovirus Receptor-Related Protein 2/PVRL2/CD112

Note: