

Human BMP Receptor II/BMPRII/PPH1



Product Name	Human BMP Receptor II/BMPRII/PPH1
Catalog No	IT-401-009p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged Human BMP Receptor II/BMPRII/PPH1 (aa 27-150; GenBank# Q13873).
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Human BMP Receptor II/BMPRII/PPH1

Note: