

Human GDNF Receptor α -2/GFRA2



Product Name	Human GDNF Receptor α -2/GFRA2
Catalog No	IT-401-012p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged Human GDNF Receptor α -2/GFRA2 (aa 22-441; GenBank# O00451).
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	$\geq 95\%$
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Human GDNF Receptor α -2/GFRA2

Note: