

# VP2 (FMDV/TypeA)



---

Product Name	VP2 (FMDV/TypeA)
Catalog No	IT-016-012Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg protein in 50mM Tris-HCl (pH8.0), 150mM NaCl and 8M urea.
Protein	6xHis tagged VP2 (FMDV/TypeA, strain A24 Cruzeiro) (aa 287-504, GenBank Accession No. P03309).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>92% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

