

# Anti-G [Hendra henipavirus]



Product Name	Anti- G [Hendra henipavirus]
Catalog No	IT-026-001M2
Source	Mouse monoclonal antibody, clone 4B6
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing glycoprotein [Hendra henipavirus] (IT-026-001p).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with glycoprotein (Hendra henipavirus)(IT-026-001p), and strongly cross with glycoprotein (Nipah henipavirus)(IT-026-002p).
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

