

Anti-nsp7 (HP-PRRSV)



Product Name	Anti- nsp7 (HP-PRRSV)
Catalog No	IT-023-002M9
Source	Mouse monoclonal antibody, clone 5D2F3
Isotype	IgG2b
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing nsp7 (HP-PRRSV) protein (IT-023-002Ep).
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with nsp7 (HP-PRRSV) as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

