


matrix protein VP40 (Mengla dianlovirus/Rousettus-wt/CHN/2015)



Product Name	matrix protein VP40 (Mengla dianlovirus/Rousettus-wt/CHN/2015)
Catalog No	IT-014-014Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg protein in 50mM Tris-HCl (pH8.0), 150mM NaCl and 8M urea.
Protein	matrix protein VP40 (Mengla dianlovirus/Rousettus-wt/CHN/2015) (aa 1-304, GenBank Accession No. AZL87825) with 6xHis tagged at C-terminus.
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>90% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: VP40 (Mengla dianlovirus/Rousettus-wt/CHN/2015)