

Human Thrombomodulin / BDCA-3 / CD141 / THBD



Product Name	Human Thrombomodulin / BDCA-3 / CD141 / THBD
Catalog No	IT-300-136p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	Human Thrombomodulin / BDCA-3 / CD141 / THBD (aa 19-515; GenBank# P07204) with 6xHis tagged at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Human Thrombomodulin / BDCA-3 / CD141 / THBD

Note: