




Human Creatine kinase B-type/CKB

Product Name	Human Creatine kinase B-type/CKB
Catalog No	IT-010-021Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg protein in 50mM Tris-HCl (pH8.0), and 150mM NaCl.
Protein	N-terminal 10xhis tagged Human Creatine kinase B-type/CKB (a.a. 1-381) (GenBank Accession No. P12277).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Human Creatine kinase B-type/CKB