


# Human Scavenger receptor cysteine-rich type 1 protein M130 / sCD163 form 4



Product Name	Human Scavenger receptor cysteine-rich type 1 protein M130 / sCD163 form 4
Catalog No	IT-300-146p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	Human Scavenger receptor cysteine-rich type 1 protein M130 / sCD163 form 4 (aa 42-369; GenBank# Q86VB7) with a 6xHis tag at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.
Reference	



E: Human sCD163 form 4

Note: