


Human ULBP2 / UL16-binding protein 2 / RAET1H / N2DL2



Product Name	Human ULBP2 / UL16-binding protein 2 / RAET1H / N2DL2
Catalog No	IT-300-150p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	Human ULBP2 / UL16-binding protein 2 / RAET1H / N2DL2 (aa 26-216; GenBank# Q9BZM5) with a 6xHis tag at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

Reference



E: Human RAET1H / N2DL2

Note: