


gp120(TRO11)(HIV/Clade B)



| | |
|----------------|--|
| Product Name | gp120(TRO11)(HIV/Clade B) |
| Catalog No | IT-001-RP4p |
| Source | Recombinant protein purified from 293 cell culture |
| Formulation | Each vial contains 100 µg of purified protein (1 mg/ml) in PBS. |
| Protein | gp120(TRO11)(HIV/Clade B) (aa 31-517; GenBank# AAW64260) with a 6xHis tag at C-terminus. |
| Applications | WB standard, antibody ELISA, Functional analysis, etc |
| Purity | ≥95% |
| Storage | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |
| Reference | |



E: gp120(TRO11)(HIV/Clade B)

Note: