


gp140(BJOX2000)(HIV/Clade CRF07)



Product Name	gp140(BJOX2000)(HIV/Clade CRF07)
Catalog No	IT-001-RP6ΔTMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	gp140(BJOX2000)(HIV/Clade CRF07) (aa 31-665; GenBank# ADI62385) with a 6xHis tag at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.
Reference	



E: gp140(BJOX2000)(HIV/Clade CRF07)

Note: