

# Anti- VP40 (Mengla dianlovirus/Rousettus-wt/CHN/2015)



Product Name	Anti- VP40 (Mengla dianlovirus/Rousettus-wt/CHN/2015)
Catalog No	IT-014-014M1
Source	Mouse monoclonal antibody, clone 4C10
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing VP40 (Mengla dianlovirus/Rousettus-wt/CHN/2015) (IT-014-014Ep).
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with VP40 (Mengla dianlovirus/Rousettus-wt/CHN/2015) as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

