


Human Cell surface glycoprotein CD200 receptor 2 / CD200R1L



Product Name	Human Cell surface glycoprotein CD200 receptor 2 / CD200R1L
Catalog No	IT-300-170p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 50 µg of purified protein (1 mg/ml) in PBS.
Protein	Human Cell surface glycoprotein CD200 receptor 2 / CD200R1L (aa 24-239; GenBank# Q6Q8B3) with 6xHis tag at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥90%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.
Reference	



E: Human Cell surface glycoprotein CD200 receptor 2 / CD200R1L

Note: