

## Pp52 (UL44) (Cytomegalovirus)



Product Name	Pp52 (UL44) (Cytomegalovirus)
Catalog No	IT-005-060Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg protein in 50mM Tris-HCl (pH8.0), 150mM NaCl and 8M urea.
Protein	The E.coli derived 51 kDa recombinant protein contains the CMV Pp52 (UL44) immunodominant regions, 202-434 amino acids (GenBank# AAS48949), fused with 26 kDa GST-tag.
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>92% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

