



## Anti- Human Creatine kinase B-type/CKB

Product Name	Anti-Human Creatine kinase B-type/CKB
Catalog No	IT-010-021M1
Source	Mouse monoclonal antibody, clone 1F10G5
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing Human Creatine kinase B-type/CKB (IT-010-021Ep).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with Human Creatine kinase B-type/CKB protein as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

