


# HA( $\Delta$ TM)(A/Wisconsin/588/2019)(H1N1)



---

Product Name	HA( $\Delta$ TM)(A/Wisconsin/588/2019)(H1N1)
Catalog No	IT-003-00117 $\Delta$ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 $\mu$ g of purified protein (1 mg/ml) in PBS.
Protein	HA( $\Delta$ TM)(A/Wisconsin/588/2019)(H1N1)(aa 18-529; GenBank# QLF80377) with 6xHis tagged at C-terminus. Most in trimer.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	$\geq$ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA( $\Delta$ TM)  
(A/Wisconsin/588/2019)(H1N1)

Note: