

Anti- Human Folate receptor alpha / FR-alpha / Ovarian tumor-associated antigen MOv18



Product Name	Anti- Human Folate receptor alpha / FR-alpha / Ovarian tumor-associated antigen MOv18
Catalog No	IT-400-029M1
Source	Mouse monoclonal antibody, clone 2D8F9
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing Human Folate receptor alpha / FR-alpha / Ovarian tumor-associated antigen MOv18.
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with Human Folate receptor alpha / FR-alpha / Ovarian tumor-associated antigen MOv18 as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

