


Spike RBD (CoVID-19/South Africa mutant)



Product Name	Spike RBD (CoVID-19/South Africa mutant)
Catalog No	IT-002-031p-SA
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 50 µg of purified protein (1 mg/ml) in PBS.
Protein	The new coronavirus Spike protein RBD (CoVID-19/South Africa mutant) (aa 339-537 having K417N, E484K, N501Y; GIAID# EPI_ISL_402121) with Human Fc and 6xHis tag at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.
Reference	



E: RBD (CoVID-19/ South Africa mutant)

Note: