

Matrix protein (Borna disease virus /strain V)



Product Name	Matrix protein (Borna disease virus/strain V)
Catalog No	IT-028-003Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg protein in 50mM Tris-HCl (pH8.0), 150mM NaCl and 8M Urea.
Protein	Matrix protein (Borna disease virus/strain V) (aa 1-142; GenBank# P0C795) with 6xHis tagged at N-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Matrix protein (BDV-V)

Note: