

Anti- NA(A/PAN)



| | |
|----------------|--|
| Product Name | Anti- NA(A/PAN) |
| Catalog No | IT-003-NA-M1 |
| Source | Mouse monoclonal antibody, clone 4D1E8 |
| Isotype | IgG2b |
| Formulation | Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol. |
| Immunogen | DNA vaccine expressing NA(A/Michigan/45/2015)(H1N1) (IT-003-1011p). |
| Applications | WB, ELISA, IP, IF etc. |
| Epitope | Linear |
| Specificity | Reacts with neuraminidase (NA) of Influenza A virus as tested; no cross with neuraminidase (NA) of Influenza B virus. |
| Storage | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |

