

Anti-Human ULBP2 / UL16-binding protein 2 / RAET1H / N2DL2



Product Name	Anti-Human ULBP2 / UL16-binding protein 2 / RAET1H / N2DL2
Catalog No	IT-300-150M3
Source	Mouse monoclonal antibody, clone 5C8G3
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing Human ULBP2 (IT-300-150p).
Applications	WB, ELISA, IP, IF, etc.
Epitope	Linear
Specificity	Reacts with Human ULBP2 / UL16-binding protein 2 / RAET1H / N2DL2 as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

