

Anti-H1(H1N1) (A/California/06/2009)



Product Name	Anti- H1(H1N1)(A/California/06/2009)
Catalog No	IT-003-001M28
Source	Mouse monoclonal antibody, clone 3A5
Isotype	IgG2a
Formulation	Each vial contains 100µg purified IgG in 100µl PBS with 40% glycerol
Immunogen	DNA vaccine expressing HA(ΔTM)(H1N1)(A/California/06/2009) protein (IT-003-SW1ΔTMp).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with HA(A/California/2009)(H1N1), not react with other strains as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

