



VP2 (Human parvovirus B19)

| | |
|----------------|--|
| Product Name | VP2 (Human parvovirus B19) |
| Catalog No | IT-015-021Ep |
| Source | Recombinant protein purified from E.coli |
| Formulation | Each vial contains 100 µg protein in 50mM Tris-HCl (pH8.0), 150mM NaCl and 8M Urea. |
| Protein | C-terminal 6xHis tagged capsid protein VP2 (Human parvovirus B19) (a.a. 1-554, GenBank# YP_004928148) |
| Applications | WB standard, antibody ELISA, Functional analysis, etc |
| Purity | ≥95% |
| Storage | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |



E: VP2 (Human parvovirus B19)

Note: