

## Anti- Human Leukemia inhibitory factor receptor / LIFR / CD118



Product Name	Anti- Human Leukemia inhibitory factor receptor / LIFR / CD118
Catalog No	IT-300-187M1
Source	Mouse monoclonal antibody, clone 2F7B8
Isotype	IgG2b
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing Human Leukemia inhibitory factor receptor / LIFR / CD118.
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with Human Leukemia inhibitory factor receptor / LIFR / CD118 as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

