

# Anti- Human MUC-1/CA 15-3/KL-6/CD227 (aa 24-240)



Product Name	Anti-Human MUC-1/CA 15-3/KL-6/CD227 (aa 24-240)
Catalog No	IT-300-189M1
Source	Mouse monoclonal antibody, clone 5C8C9
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing Human MUC-1/CA 15-3/KL-6/CD227 (GenBank P15941-1; aa 24-240).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with Human MUC-1/CA 15-3/KL-6/CD227 (aa 24-240) protein as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

