


Spike RBD (CoVID-19/BA.2.75 variant), His-tagged



Product Name	Spike RBD (CoVID-19/BA.2.75 variant), His-tagged
Catalog No	IT-002-036p-BA.2.75
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 50 µg of purified protein (1 mg/ml) in PBS.
Protein	The new coronavirus Spike protein RBD (CoVID-19/BA.2.75 variant) (G339H, S371F, S373P, S375F, T376A, D405N, R408S, K417N N440K G446S N460K S477N T478K E484A Q498R N501Y Y505H; GIAID# EPI_ISL_402121) with 6xHis tag at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.
Reference	



E: RBD (CoVID-19/BA.2.75), his-tagged

Note: