


# EEV glycoprotein A35R [Monkeypox virus]



|                |  |
|----------------|--|
| Product Name   | EEV glycoprotein A35R [Monkeypox virus]  |
| Catalog No     | IT-012-015p  |
| Source         | Recombinant protein purified from 293 cell culture   |
| Formulation    | Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.  |
| Protein        | EEV glycoprotein A35R [Monkeypox virus] (aa 57-181, Genbank No. YP_010377142) with a 6xHis tag at C-terminus.  |
| Applications   | WB standard, antibody ELISA, etc   |
| Purity         | ≥ 95%  |
| Storage        | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures.   |



E: EEV glycoprotein A35R [Monkeypox virus]

Note: