


IMV membrane protein L1R [Monkeypox virus]



Product Name	IMV membrane protein L1R [Monkeypox virus]
Catalog No	IT-012-016p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	IMV membrane protein L1R [Monkeypox virus] (aa 2-183, Genbank No. YP_010377077) with a 6xHis tag at C-terminus.
Applications	WB standard, antibody ELISA, etc
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: IMV membrane protein L1R [Monkeypox virus]

Note: