



gp120(BL035.W6M.ENV.C1)

Product Name	gp120(BL035.W6M.ENV.C1)(Clade A)
Catalog No	IT-001-115p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HIV-1 gp120 (BL035.W6M.ENV.C1)(Clade A) protein (a.a.35-499) (Genbank No.DQ208480).
Applications	WB standard, antibody ELISA, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

