

Anti-H1(A/Christchurch/16/2010)(H1N1)



Product Name	Anti-H1(A/Christchurch/16/2010)(H1N1)
Catalog No	IT-003-00103M1
Source	Mouse monoclonal antibody, clone ITF7
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	Recombinant HA(ΔTM)(A/Christchurch/16/2010)(H1N1) protein (IT-003-00103ΔTMp).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with 2013-2014 seasonal flu H1(A/Christchurch/16/2010) (H1N1). Not tested for its cross-reactivity to other flu subtypes.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-003-00103M1

