

Anti-HA(H1N1/2009 Swine Flu)



Product Name	Anti-HA(H1N1/2009/Swine Flu)
Catalog No	IT-003-001M6
Source	Mouse monoclonal antibody, clone 39F1
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS.
Immunogen	Recombinant HA(A/California/06/2009(H1N1)) protein (a.a. 18-530) (Genebank accession ACP41935).
Applications	ELISA, WB, IP, IF etc.
Epitope	Linear
Specificity	Reacts with HA(H1N1/2009/Swine Flu) only. Minimal cross-reactivity to other subtypes.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

