

# Anti-HA1(A/Texas/50/2012)(H3N2)



Product Name	Anti-HA1(A/Texas/50/2012)(H3N2)
Catalog No	IT-003-00427M1
Source	Mouse monoclonal antibody, clone ITE9
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	Recombinant HA1 proteins (H3N2)(A/Texas/50/2012)(IT-003-00427p)
Applications	ELISA, WB, IP, IF etc.
Epitope	Linear
Specificity	Reacts with 2013-2014 seasonal flu HA1(A/Texas/50/2012)(H3N2). Not tested for its cross-reactivity to other subtypes.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

**ELISA Profile of IT-003-00427M1**

