




HA1 (H₅N₁)(A/Hong Kong/213/2003)

Product Name	HA1 (H ₅ N ₁)(A/Hong Kong/213/2003)
Catalog No	IT-003-0059p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA1(H ₅ N ₁)(A/Hong Kong/213/2003) protein (a.a.17-539; Genbank# BAE07201).
Applications	Western Blotting standard, antibody ELISA, immunogen, etc.
Purity	≥ 95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA1 (H₅N₁)(A/Hong Kong/213/2003)