

Anti-N2(H3N2)



Product Name	Anti-N2 (H3N2)
Catalog No	IT-003-020M3
Source	Mouse monoclonal antibody, clone 218C11
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	Recombinant NA protein (A/Wisconsin/67/X-161/2005(H3N2)) (Genebank accession No. ABO37609)
Applications	WB
Epitope	Linear
Specificity	Reacts with most seasonal strain N2 (H3N2). Cross-reacts with N2 (H9N2). Cross-reactivity to other subtypes not tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

