

HA1(A/flat-faced bat/Peru/033/2010)(H18N11)



Product Name	HA1(A/flat-faced bat/Peru/033/2010)(H18N11)
Catalog No	IT-003-050p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA1(A/flat-faced bat/Peru/033/2010)(H18N11) protein (a.a.15-339; Genbank# CY125945).
Applications	Western Blotting standard, antibody ELISA, immunogen, etc.
Purity	≥ 95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

