


# NS1 (DV3)



---

Product Name	NS1 (DV3)
Catalog No	IT-006-0052p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (100 µg/350 µl) in PBS.
Protein	C-terminal 6xHis tagged Dengue Virus-3 NS1 protein + V5 (Genebank No. FJ744738.1).
Applications	WB standard, antibody ELISA, vaccine, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: NS1 (DV3) protein

Note: