

Anti-NS1(Flavivirus)



Product Name	Anti-NS1(Flavivirus)
Catalog No	IT-006-NS1-Pan-M1
Source	Mouse monoclonal antibody, clone 1F21
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	DNA vaccine expressing NS1(Tick-borne encephalitis virus) (a.a. 1-371) (Genebank No. NP_043135)
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with NS1 from all Flavivirus members as tested.
Storage	Store at 4°C, stable for 6months from the date of shipment. For long term storage, aliquot, and freeze at -80°C. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

