




GP (Ebola/Sudan - Nakisamata)

Product Name	GP (Ebola/Sudan - Nakisamata)
Catalog No	IT-014-007p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	Ebola [Sudan – Nakisamata] glycoprotein (aa 1-647) (Genbank No. AFP28231)
Applications	WB standard, antibody ELISA, etc
Purity	≥ 95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: GP1/GP2 (Ebola/Sudan - Nakisamata) proteins

Note: