


diphtheria toxin [*Corynebacterium diphtheriae*]



Product Name	diphtheria toxin [<i>Corynebacterium diphtheriae</i>]
Catalog No	IT-505-001Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg protein in 50mM Tris-HCl (pH8.0), 150mM NaCl and 8M Urea.
Protein	C-terminal 6xHis tagged diphtheria toxin [<i>Corynebacterium diphtheriae</i>] (a.a. 26-560, GenBank#WP_072564851)
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥90%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: diphtheria toxin [*Corynebacterium diphtheriae*]

Note: