


Capsid protein VP2 (Tremovirus A/AEV Calnek)



Product Name	Capsid protein VP2 (Tremovirus A/AEV Calnek)
Catalog No	IT-016-032Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg protein in 50mM Tris-HCl (pH8.0), 150mM NaCl and 8M Urea.
Protein	C-terminal 6xHis tagged Capsid protein VP2 (Tremovirus A/AEV Calnek) (a.a. 20-242, GenBank#Q9YLS4)
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: VP2 (Tremovirus A/AEV Calnek)

Note: