


gp140(O_MVP5180)(HIV/Clade O)



Product Name	gp140 (O_MVP5180) (HIV/Clade O)
Catalog No	IT-001-196p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	gp140 (HIV-1/O_MVP5180) (aa 34-693; GenBank# Q79670) with a 6xHis tag at C-terminus. Most are in trimer form
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥90%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.
Reference	



E: gp140 (HIV-1/O_MVP5180)

Note: