

# Anti-EEV glycoprotein A35R [Monkeypox virus MPXV]



Product Name	Anti-EEV glycoprotein A35R [Monkeypox virus MPXV]
Catalog No	IT-012-015M2
Source	Mouse monoclonal antibody, clone 4E10D11
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing EEV glycoprotein A35R [Monkeypox virus MPXV] (IT-012-015p).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with EEV glycoprotein A35R [Monkeypox virus MPXV] as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

