

## Anti-NP(PAN)



Product Name	Anti-NP(PAN)
Catalog No	IT-003-048M7
Source	Mouse monoclonal antibody, clone 2D9C5
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing NP(B/Wisconsin/01/2010) (IT-003-048Ep).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with NP(B/Wisconsin/01/2010) and NP(A/PR/8/1934) (H1N1)) strongly, with NP(A/Brisbane/10/2007 (H3N2)) weakly.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

### 2D9C5

