

Anti-IEEV membrane glycoprotein B6R [Monkeypox virus]



Product Name	Anti-IEEV membrane glycoprotein B6R [Monkeypox virus]
Catalog No	IT-012-017M8
Source	Mouse monoclonal antibody, clone 4A6E5
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing IEEV membrane glycoprotein B6R [Monkeypox virus] (IT-012-017p).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with IEEV membrane glycoprotein B6R [Monkeypox virus] as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

