

capsid protein (Norovirus Hu/GII-4/CGMH23/2010/TW)



Product Name	capsid protein (Norovirus Hu/GII-4/CGMH23/2010/TW)
Catalog No	IT-016-021p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	capsid protein (Norovirus Hu/GII-4/CGMH23/2010/TW) (a.a. 267-418) (GenBank Accession No. AFN09984) with 6xHis tag at C-terminal.
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

