

gp120(NL4-3)(HIV-1/Clade B)



Product Name	gp120(NL4-3)(HIV/Clade B)
Catalog No	IT-001-204p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	gp120(NL4-3)(HIV/Clade B) (aa 33-505; GenBank#Q72502) with a 6xHis tag at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

Reference



E: gp120(NL4-3)(HIV/Clade B)

Note: